VeriPRO Quick Reference Guide

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**HOWARD
LEIGHT** 
by Honeywell

Welcome

Howard Leight's VeriPRO makes it easy to get an accurate, real-world picture of your employees' hearing protection.

This Quick Reference Guide gives you an overview of VeriPRO hardware and software and gets you on your way to providing a valuable hearing conservation service to your employees. The guide provides information needed to set up, install and begin using VeriPRO.

1.0 VeriPRO at a Glance

Developed by Howard Leight, in conjunction with the House Ear Institute, VeriPRO makes it easy to get an accurate, real-world picture of your employees' hearing protection. Find out whether your employees are receiving optimal protection, require additional training on how to fit their earplugs, or need to try a different model. VeriPRO uses sophisticated software in a user-friendly format to measure the Personal Attenuation Rating [PAR] – the real-ear attenuation your employees are receiving from their earplugs.

Designed to be as easy to use as an automated cash machine or self-service ticket kiosk, VeriPRO is a simple yet powerful component of any Hearing Conservation program.

- Measures the Personal Attenuation Rating (PAR) for each employee
- Captures real-world attenuation by using unmodified earplugs
- Works with any earplug, corded or uncorded
- Fulfills regulatory requirements to ensure proper initial fitting of hearing protection
- Creates accurate, easy-to-understand reports within minutes
- Provides an ideal opportunity for one-on-one education


1.1 How It Works

Unlike typical audiometric tests that measure the softest level at which a listener can hear a tone, VeriPRO utilizes a relative measurement of hearing level – the loudness of a tone in one ear is balanced to a tone in the opposite ear at a known level. This unique “loudness-balance” method allows an accurate hearing test to be performed with normal background noise levels found in a typical office or worksite – no special sound booth or accommodations are needed.

The VeriPRO test procedure

VeriPRO’s test method is simple to learn and fast to complete. Users balance tones using an on-screen slider bar so that the tones are equally loud in both ears.

1. VeriPRO measures the user’s relative hearing level without earplugs, across a range of frequencies.
2. The user inserts an earplug in the right ear and repeats the loudness-balance procedure.
3. The user repeats the loudness-balance procedure with earplugs in both ears.

1.1 How It Works (continued)

Personal Attenuation Rating (PAR)

Based on the results, VeriPRO calculates the amount of attenuation at each frequency in each ear and displays a Personal Attenuation Rating (PAR) for each ear. The users know exactly how much attenuation they are getting from their earplugs the instant they are fitted *(see Section 3.3)*.

Quick Check Results

Quick Check uses the same loudness balance procedure, but only tests one critical frequency (500 Hz) in each ear *(see Section 3.4)*.

1.2 VeriPRO Components

CD

Includes software, specially designed interface for users and reporting features.

Audio Processor

Powers the system and serves as an external sound card and headphone amplifier. Plug the audio processor into the USB port of your computer and plug the headphones into the audio processor.

Headphones

Audiometrically balanced specifically for use with VeriPRO, these headphones provide optimal audio output and greater backset of the speaker for testing earplugs. Color-coded earcups: Red for right, Blue for left.

Note: *Do not use other headphones with VeriPRO, as this will affect test results and may cause discomfort if the earplugs make contact with the headphone.*

1.3 System Requirements

Minimum Requirements

700Mhz CPU, 256MB RAM

Recommended Requirements

1.2Ghz CPU or faster,
512MB RAM

Hard-Disk Space Requirements

500MB

Software Requirements

Windows XP, Windows 7 or
Windows Vista with SP3 or later,
Microsoft.NET Framework 2.0

When using Windows 7 or Windows Vista, VeriPRO must be installed/opened with Administrator privileges.

2.0 Getting Started

To get started using VeriPRO, simply plug in the hardware components, load the software onto your computer and enter users into the database.

1. Plug VeriPRO audio processor into USB port of your computer
2. Plug VeriPRO headphones into audio processor
3. Turn off computer-generated sounds in your Control Panel
4. Install VeriPRO software
5. Enter user names into the database
6. Run VeriPRO test

2.1 Install & Register Software

To install the VeriPRO software onto your computer:

1. Insert the **VeriPRO** CD into your computer's CD drive.
2. Double-click on the **My Computer** icon on your desktop or in your Start menu.
3. Access your CD drive by clicking on the D:\ drive icon.
4. Double-click on the **VeriPRO** icon to begin the installation process.
5. Follow the on-screen instructions to install VeriPRO in the directory of your choice.

2.2 Managing Users

You must add user names and other information into the database before you begin testing.

To add users to the database:


1. Open program by clicking the **VeriPRO** icon.
2. At the welcome screen, click Continue. This brings you to the Main Menu screen.
3. Click the **Toolbox** icon in the top right corner.
4. Enter the **Administrator Password** into the pop-up box. Type the word **"password"** for your initial log-in, then change the password to your own preference. This password is encoded and cannot be retrieved by anyone but you.
5. Click **Manage Users**.

2.2 Managing Users (continued)

To enter user names:

1. Click **Add User**.
2. Enter an employee's first / middle / last name, ID number, Noise Exposure Level, Department, and Work Location. Required fields are indicated with an asterisk (*). After a Department or Work Location has been entered initially, these entries can be selected from the drop-down menu in future tests. They do not need to be entered again. Instead of entering user names individually, you can also upload a batch file of users. Click on Upload Users and follow the instructions to prepare a tab-delimited text file, with data in appropriate columns. Once the file is prepared, click on Upload File.
3. Enter the employee's noise exposure level, if known. This exposure level should be the 8-hour Time-Weighted Average exposure level, and may be entered as dBA (most common) or dBC. This entry must be between 70 - 170 dB.

2.3 Set-Up and Preferences

Preferences

This function allows the administrator to select various testing options.

1. Select your preferred language.
2. Select the rating method used in your region (NRR for North America, SNR for Europe, SLC₈₀ for Australia).
3. Select your Safe Noise Criterion. This is the level you would like to define as a safe noise exposure. We recommend 80 dB as the Safe Noise Criterion; however, some companies and regulatory agencies use other levels between 80-85. Your selection of Safe Noise Criterion does not affect the VeriPRO test itself, but does affect the calculation of Safe Exposure Level in the results.


The other settings in the Preferences screen allow the administrator to adjust the parameters of the VeriPRO test, including audio output levels and timing. (continued on next page)

VeriPRO

HOWARD
LEIGHT
by Honeywell

2.3 Set-Up and Preferences (continued)

Typically no adjustment is needed in these settings – the default settings accommodate nearly all users. The settings are programmed as shown:

Upload Company Logo

VeriPRO allows you to customize the screen with your company logo. Simply upload a .jpg, .gif or .png file. Image dimensions should not exceed 416 pixels W x 33 pixels H.

Quick Reference Guide

An electronic copy of this Quick Reference Guide is available from this menu.

Export Data

VeriPRO results may be exported to an external file for data collection and analysis.

1. In the **Toolbox** menu, select **Export Data**.
2. Select the date range for the tests you would like to export.
3. Click on **Export**.

2.4 Audio Settings

Before you begin running tests, you must disable the internal computer-generated sounds in the Windows Sound Control Panel. Otherwise, computer sounds (such as incoming e-mail alerts) will sound uncomfortably loud for users and will interfere with the test.

Disable internal sounds

From your Windows desktop, click:


- **Start**
- **Settings**
- **Control Panel**
- **Sounds and Audio Devices**

2.4 Audio Settings (continued)

From the Sounds and Audio Devices window

- Open the **Sounds** tab
- From the Sound Scheme window, select **No Sounds** and then click "OK"

Note: Simply turning down the volume control or muting the volume from the Windows desktop does NOT turn off computer-generated sounds.

3.0 Using VeriPRO

VeriPRO is designed to be easy to use, requiring minimal training for you and your employees. It provides options for Complete Check and Quick Check, as well as a variety of reporting formats for administrators.

3.1 Main Menu

Complete Check

Leads users through a sequence of five test frequencies in three different parts. End result is a Personal Attenuation Rating (PAR). Usually completed in eight to ten minutes.

Quick Check

Ideal for earplug selection, training and periodic testing. Tests at one critical frequency: 500 Hz. Usually completed in three minutes or less.

Reports

Provides options for viewing and printing VeriPRO test results. Reports can include a selected user's most recent individual test, any past individual test, or a summary of all VeriPRO tests on file for a given user.

Fit Training

Helpful on-screen instructions and videos demonstrating proper fit of Howard Leight earplugs.

3.2 Select Earplug Screen

From the Main Menu, users are directed to the Earplug Selection screen. Although Howard Leight® earplugs are shown in the selection screen, a user can perform a VeriPRO test with other types of earplugs by selecting **Other Earplugs**.

If using 'Other Earplug', enter the model name and rated attenuation value of that earplug. After an earplug has been entered the first time, it can be selected from the drop-down menu.

3.3 Complete Check

Complete Check leads the user through a sequence using five test frequencies in three different parts:

Part 1 – Loudness balance with *no earplugs* (Both Ears Unoccluded)

Part 2 – Loudness balance with earplug in *right ear only* (Right Ear Occluded) **NOTE:** After Part 2 is complete, it is important to keep the right earplug in place. Remind users to not remove the right earplug after Part 2.

Part 3 – Loudness balance with earplugs in *both ears* (Both Ears Occluded)

The frequencies used for each test are 500 Hz (repeated twice as a reliability check), 250 Hz, 1000 Hz, 2000 Hz and 4000 Hz. Complete Check is ideal for use as a periodic check for workers with a significant shift in their hearing.

Quality Controls

To increase test accuracy, four quality checks are built into the VeriPRO software:

Roving Slider Bar – An arbitrary scale that changes with every new screen. This prevents users from “benchmarking” their loudness-balance, as the middle position of one slider bar will sound very different from the middle position of the next slider bar.

500 Hz Retest – Demonstrates retest reliability within 6 dB between the first and second trials.

Minimum Attenuation – Warns users if they are not achieving a minimum level of attenuation. Users can refit the earplug if they wish. If minimum attenuation is still not achieved, the user will be warned a second time and the test will continue.

Slider Bar Movement – Tracks how many times the slider bar is engaged. If a user leaves the slider bar in the same position throughout the test, a warning is displayed on the screen and the Results screen encourages the user to retest.

3.3.1 Using the Slider Bar

To proceed with the VeriPRO test, the user moves the slider bar up and down to adjust the loudness of the left ear until it matches the reference tone in the right ear. When volume in the left ear and right ear match, press MATCH.

- To extend the range of the slider bar higher or lower, press the double-arrows at the top or bottom of the slider bar. Note that the double arrows flash red when the slider bar range has been extended.
- If user cannot match the loudness of the tones, press CANNOT MATCH and click on the appropriate remark.

Remark	Effect
I couldn't hear the sound in my left ear.	Clicking on this choice adjusts the dB range 10 dB higher in the left ear.
I couldn't hear the sound in my right ear.	Clicking on this choice ends the test, since the subject cannot hear even the reference tone. To correct this, the System Operator should adjust the Reference Ear Level in the Both Ears Unoccluded section of Preferences.
I couldn't hear the sound in either ear.	The test ends, similar to above.
Sound in my left ear is softer than in my right ear.	Clicking on this choice re-adjusts the scale on the slider bar 10 dB higher.
Sound in my left ear is louder than in my right ear.	Clicking on this choice re-adjusts the scale on the slider bar 10 dB lower.

3.4 Quick Check



The Quick Check mode is a fast screening check of the earplug fit. It is ideal for earplug selection, training and periodic testing. It uses the same test procedures as the Complete Check mode, but only tests at one critical frequency: 500 Hz. If an earplug is fit poorly enough to allow noise to easily leak through, this always shows up at 500 Hz.

The Quick Check leads the user through the same three parts: loudness balance with no earplugs, with earplug in right ear only, and with earplugs in both ears. Due to testing at only one frequency instead of five frequencies, the test is completed much more quickly – usually in three minutes or less for most users – but is not as accurate as the Complete Check. The result of the Quick Check is simply the attenuation of the earplug at 500 Hz.

3.5 User Reports


Under this menu, you can select the way you would like to view VeriPRO reports:

- Select the test you would like to view
- Select optional views of the report (attenuation chart by frequency or summary of all tests)
- Select whether you would like to create a PDF file of the report to save on your computer, or view as a printable document.

3.6 Earplug Fit Training

VeriPRO offers helpful on-screen fitting instructions and fitting videos for all Howard Leight earplugs. Users select which style they are using. The fit training video lasts less than two minutes and provides step-by-step text instructions and video demonstrations on how to insert an earplug properly for optimal attenuation.

4.0 Troubleshooting

VeriPRO System	Problem	Recommendation
Overall	Problems with hardware, software loading, computer compatibility.	Contact VeriPRO Technical Support.
Headphones	Hissing or popping sound in headphone.	Due to the high output levels required for the VeriPRO headphones, some background hissing or popping sounds are normal. Since VeriPRO tests above threshold levels, these headphone noises do not affect the test results.
Headphones	The tones are heard in both ears at the same time, instead of just one ear at a time.	Reconnect the headphone cable to the audio processor. Make sure the cable is fully inserted — partial insertion will cause the signal to be heard in both channels simultaneously.
Audio Processor	Power LED display does not light.	The LED will light whenever power is going into the audio processor. If the LED is not on, try reconnecting the mini-USB cable to the USB port on your computer and to the audio processor.
Software	When program launches, there is a warning box: System May Not Be Properly Calibrated.	VeriPRO cannot identify the correct audio processor or headphone. <ul style="list-style-type: none">• Make sure only the VeriPRO audio processor and VeriPRO headphones are used with the software.• Check all cable connections. Reconnect the headphone cable to the audio processor. Make sure the cable is fully inserted — partial insertion will cause the signal to be heard in both channels simultaneously.
Sounds	No tones are audible.	<ul style="list-style-type: none">• Check all cable connections.• Restart the VeriPRO software, allowing your computer to recognize the VeriPRO hardware ["New Hardware"] to be found by the operating system.

4.0 Troubleshooting (continued)

VeriPRO System

Problem

Recommendation

Sounds

The tones are too loud.

There are four reasons why tones in the headphone may sound too loud:

1. System Operator did not disable internally-generated computer sounds (see Section 2.4). Make sure the Sound Scheme is set to "No Sounds" during VeriPRO testing.
2. Headphone was switched. Make certain the red earcup (R) is on the right ear and the blue earcup (L) is on the left ear.
3. Earplugs were not inserted. The VeriPRO software generates loud tones during those parts of the test for which earplugs are worn. Make sure earplugs are inserted correctly at appropriate times during the test sequence.
4. Other programs are running on the computer's desktop. For best results, do not launch or run other programs while testing with VeriPRO. These can interfere with the test and the presentation of tones.

Sounds

The background noise in the test room is so loud, users cannot hear the VeriPRO tones.

VeriPRO can compensate for loud background noise:

- Open the System Toolbox and click on Preferences to access the settings.
- Raise the Reference Ear Level (Both Ears Unoccluded) to 70 dB.
- The reference tone will be easier to hear over the background noise without affecting validity of test results.

4.1 Care and Maintenance

Daily Listening Check

To check headphone output, the VeriPRO administrator should perform a short listening check prior to testing:

- Connect all hardware and open the VeriPRO program.
- Start a Complete Check at 500 Hz and balance the loudness of the tones. For a listener with normal hearing, if several Cannot Match responses are required before you reach a loudness balance, the headphone output may be questionable. A calibration check is recommended.

Cleaning of Headphones

Recommendations for cleaning the VeriPRO headphones are similar to other earmuffs. Wipe down the cushions and headband with mild soap and damp cloth, and dry before use. Do not submerge headphones in water. Replace the cushions if they appear stiff or cracked. Extra cushions are included with each VeriPRO system, and can also be ordered from your safety distributor.

Calibration Check

We recommend an annual calibration check of the VeriPRO audio processor and headphones. Contact your Howard Leight representative for details.

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About Us

Howard Leight is one of the largest global manufacturers of hearing protection in the industrial market and the recognized innovator in protection and people-oriented fit. Our products and technologies empower people to ensure the right fit and protection for every user, every environment. Visit us online at www.howardleight.com.

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www.howardleight.com
www.howardleightveripro.com

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Honeywell Safety Products USA, Inc.